TESLA SMART THERMOMETER

USER MANUAL
Model: UFR102

EN

DE

CZ

SK

PL

HU

SL

Introduction

- Your new Infrared thermometer uses advanced infrared (IR) technology to measure temperature instantly and accurately on the body or object.

- Easy to use and less measurement time

Tesla Smart Thermometer does not need to contact body or object to ensure the safety and hygiene. Its ergonomic design makes this thermometer be simple and very easy to check the temperature. It only takes 1 second to take measurement and reading.

- Body, ear and object temperature modes

Tesla Smart Thermometer supports to measure temperature of body surface, ear canal and object. Measuring range of object-temperature mode is 32.0°F-199.9°F (0.0°C -100.0°C). This means except body surface, ear canal temperature, it also

supports to take the surface temperature of object temperature as below,

- Surface temperature of milk in the bottle
- Surface temperature of baby bath
- Environment temperature
- Object mode is for reference only

- Alarm sounds

When body temperature is over 99.5°F(37.5°C), an alarm sounds for 10 times

- Memory function

It stores up to 20 sets recent measurement data.

- Auto power off

Automatically power off if left idle for 60±5seconds.

- Measurement

Measurement time 1 second and the measurement distance is within 1 cm.

Please read the manual carefully before you use the unit and keep for future reference.

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