BRAUN

ThermoScan®

Ear thermometer

IRT 6020
IRT 6520

| | | |
|-----------|-------------------|-----|
| GB | English | 3 |
| CZ | Česky | 14 |
| DK | Dansk | 25 |
| EE | Eesti | 36 |
| FI | Suomi | 47 |
| HR | Hrvatski | 58 |
| HU | Magyar | 69 |
| LT | Lietuvių | 80 |
| LV | Latviešu | 91 |
| NO | Norsk | 102 |
| PL | Polski | 113 |
| RU | Русский | 124 |
| SE | Svenska | 135 |
| SI | Slovenščina | 146 |
| SK | Slovenčina | 157 |

This product is manufactured by Kaz Europe Sàrl under a license to the 'Braun' trademark. 'Braun' is a registered trademark of Braun GmbH, Kronberg, Germany. ThermoScan is a registered trademark of Kaz Europe Sàrl.

Kaz Europe Sàrl
Place Chauderon 18
CH-1003 Lausanne
Switzerland

© 2014, Kaz Europe Sàrl

www.hot-europe.com


Patents: www.kaz.com/patents/braun
Made in Mexico.



Printed in USA

IRT6020/6520 NOEE
P/N: 31IM65NO190
04MAR14

IRT 6520

IRT 6020

Product description

IRT 6520

- 1 Lens filter (Box of 20)
- 2 Probe tip
- 3 Probe
- 4 ExacTemp light
- 5 Lens filter ejector
- 6 Display
- 7 Age Precision and Memory buttons
- 8 Power button
- 9 Start button
- 10 Battery door
- 11 Protective case

IRT 6020

- 1 Lens filter (Box of 20)
- 2 Probe tip
- 3 Probe
- 4 ExacTemp light
- 5 Lens filter ejector
- 6 Display
- 7 Power button
- 8 Start button
- 9 Battery door
- 10 Protective cap

The Braun ThermoScan thermometer has been carefully developed for accurate, safe and fast temperature measurements in the ear.

The shape of the thermometer probe prevents it from being inserted too far into the ear canal, which can hurt the eardrum.

However, as with any thermometer, proper technique is critical to obtaining accurate temperatures. Therefore, read the instructions carefully and thoroughly.

The Braun ThermoScan thermometer is indicated for intermittent measurement and monitoring of human body temperature for people of all ages. It is intended for household use only.

Use of IRT 6520 along with Age Precision™ feature is not intended as a substitute for consultation with your physician.

WARNINGS AND PRECAUTIONS

- Please consult your doctor if you see symptoms such as unexplained irritability, vomiting, diarrhea, dehydration, changes in appetite or activity, seizure, muscle pain, shivering, stiff neck, pain when urinating, etc. in spite of any color of the background light or absence of fever.
- Please consult your doctor if the thermometer shows elevated temperature (yellow or red background lights for IRT 6520).

The operating ambient temperature range for this thermometer is 10 – 40 °C (50 – 104 °F). Do not expose the thermometer to temperature extremes (below -25 °C / -13 °F or over 55 °C / 131 °F) or excessive humidity (>95% RH). This thermometer must only be used with genuine Braun ThermoScan lens filters (LF 40).

To avoid inaccurate measurements always use this thermometer with a new, clean lens filter attached.

If the thermometer is accidentally used without a lens filter attached, clean the lens (see «Care and cleaning» section). Keep lens filters out of reach of children.

This thermometer is intended for household use only. This product is not intended to diagnose any disease, but is a useful screening tool for temperature. Use of this thermometer is not intended as a substitute for consultation with your physician.

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