

BRAUN

No touch + touch Thermometer




BNT300

GB	English	2
CZ	Česky.....	11
DK	Dansk	20
FI	Suomi	29
NO	Norsk.....	38
PL	Polski.....	47
SA	اللغة العربية	64
SE	Svenska	65
SK	Slovenčina.....	74

Product description

- 1 Power button
- 2 Scanner
- 3 LCD display
- 4 Temperature button
- 5 Battery cover
- 6 Protective scanner cap
- 7 Silent mode switch
- 8 Guidance light



Certain trademarks used under license from
The Procter & Gamble Company or its affiliates.

Patents: https://www.braunhealthcare.com/us_en/patents

 Kaz Europe Sàrl
A Helen of Troy Company
Place Chauderon 18
CH-1003 Lausanne
Switzerland



KAZ Hausgeräte GmbH
Kölner Strasse 23
DE-57439 Attendorn - Germany

HOT (UK) LTD
PO Box 5984
Sheffield, S2 9GN - United Kingdom

© 2021, All rights reserved

www.braunhealthcare.com/us_en/patents



REF BNT300EE
P/N: A004442R3
02NOV21

Thank you for purchasing the Braun No touch + touch thermometer (BNT300). This thermometer is a high-quality product incorporating the latest technology and tested in accordance with international standards. With its unique technology, the product can provide a stable, accurate reading with each measurement. The thermometer performs a self-test every time it is switched on to always guarantee the accuracy of measurements.

Please read these instructions carefully before using this product and keep the instructions and the thermometer in a safe place.

Intended use and contraindications

This No touch + touch thermometer is an handheld, battery-powered, infrared thermometer intended for the intermittent measurement of human body temperature in a home-use environment on people of all ages, except pre-term and small-for-gestational-age babies, using the center of the forehead as the measurement site.

These thermometers are not intended for clinical use in a professional environment and are for home usage only. This thermometer is not intended to interpret hypothermic temperatures.



WARNINGS AND PRECAUTIONS

Keep out of reach of children under 12 years.

Never use the thermometer for purposes other than those it has been intended for. Please follow the general safety precautions when using on children.

Never immerse the thermometer into water or other liquids (not waterproof). For cleaning and disinfecting please follow the instructions in the "Care and cleaning" section.

Do not store this thermometer in temperature extremes below -25 °C or over 55 °C (below -13 °F or over 131 °F) or in excessive humidity (above 95% non-condensing relative humidity).

If thermometer is stored in a location that is cooler or warmer than where it is being used, let it sit in the patient's room for 10 minutes before taking the measurement.

Do not use the thermometer if there are signs of damage on the scanner or on the thermometer itself. If damaged, do not attempt to repair the product.

Never scratch the sensor with a sharp object when cleaning.

This thermometer consists of high-quality precision parts. Do not drop the instrument. Protect it from severe impact and shock. Do not twist the instrument or the measuring sensor.

This thermometer is intended for household use only.

Use of this thermometer is not intended as a substitute for consultation with your doctor.

Temperature elevation may signal a serious illness, especially in neonates and infants, or in adults who are old, frail, or have a weakened immune system. Please seek professional advice immediately when there is a temperature elevation and if you are taking temperature on:

- Neonates and infants under 3 months (Consult your doctor immediately if the temperature exceeds 37.4 °C [99.4 °F])
- Patients over 60 years of age (Fever may be blunted or absent in older patients)
- Patients having diabetes mellitus or a weakened immune system (e.g., HIV positive, cancer chemotherapy, chronic steroid treatment, splenectomy)

- Patients who are bedridden (e.g., nursing home patient, stroke, chronic illness, recovering from surgery)
- A transplant patient (e.g., liver, heart, lung, kidney).

Do not allow children to take their temperatures unattended.

Please consult your doctor if you see symptoms such as unexplained irritability, vomiting, diarrhea, dehydration, changes in appetite or activity, seizure, muscle pain, shivering, stiff neck, pain when urinating, etc., even in the absence of fever.

Even in the absence of fever, those who exhibit a normal temperature may still need to receive medical attention. People who are on antibiotics, analgesics, or antipyretics should not be assessed solely on temperature readings to determine the severity of their illness.

Do not modify this equipment without the authorization of the manufacturer.

Why Braun No touch + touch thermometer?

Measurement in under 2 seconds

The innovative infrared technology allows measurement of touch temperature in under 2 seconds whether you use the touch or no touch option.



Accurate and reliable


Due to the unique patented technology, the No touch + touch thermometer captures the heat naturally given off by the forehead to calculate body temperature value. Get the same professional accuracy whether touching your child's forehead or holding up to 2.5 cm away.



Easy to use

The No touch + touch thermometer is non-invasive. A measurement can be taken even while a child is sleeping.

The No touch + touch thermometer is less threatening to a child than a rectal thermometer and easier to use than other methods.



Safe and hygienic

No-touch option helps minimize spreading of germs. Safe for use on children and adults.

How does Braun No touch + touch work?

The Braun No touch + touch thermometer measures infrared energy radiated from the skin at the center of the forehead area. This captured energy is collected through the sensor and converted to an oral equivalent temperature value.

The No touch + touch thermometer has been clinically tested and proven to be safe and accurate when used in accordance with its operating instruction manual.

Obsah je uzamčen

Dokončete, prosím, proces objednávky.

Následně budete mít přístup k celému dokumentu.



Proč je dokument uzamčen? Nahněvat Vás rozhodně nechceme. Jsou k tomu dva hlavní důvody:

- 1) Vytvořit a udržovat obsáhlou databázi návodů stojí nejen spoustu úsilí a času, ale i finanční prostředky. Dělali byste to Vy zadarmo? Ne*. Zakoupením této služby obdržíte úplný návod a podpoříte provoz a rozvoj našich stránek. Třeba se Vám to bude ještě někdy hodit.

**) Možná zpočátku ano. Ale vězte, že dotovat to dlouhodobě nelze. A rozhodně na tom nezbohatneme.*

- 2) Pak jsou tady „roboti“, kteří se přiživují na naší práci a „vysávají“ výsledky našeho úsilí pro svůj prospěch. Tímto krokem se jim to snažíme překazit.

A pokud nemáte zájem, respektujeme to. Urgujte svého prodejce. A když neuspějete, rádi Vás uvidíme!