BabyScale B5 BT

EN - MANUAL

DE - ANLEITUNG

CS - NÁVOD

SK - NÁVOD

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SL - NAVODILA

HR - UPUTE

EN – MANUAL

INTRODUCTORY INFORMATION

- Before switching on the scale for the first time, insert the 3 included AAA batteries into the compartment at the bottom of the scale. Make sure they are oriented according to the diagram in the compartment.
- Place the scale on a flat, solid surface or a mat in a dry environment.
- If you use the scale on a soft surface or carpet, the product will not be able to measure the child's weight optimally.
- Do not exceed the maximum weight (20 kg).
- When not in use, do not place anything on top of it or put unnecessary strain on it.
- Do not shake the scale and avoid impacts.
- Avoid contact with water.
- For the rest of the information regarding the safe use of the device, see the Safety Instructions in the box.

FIRST USE, WEIGHING WITHOUT THE APP

1. Press the **ON/OFF/TARE** button, the scale will turn on, and the display will light up and show „Hi“.
2. When the text disappears and the scale display shows 0, the scale is ready for use.
3. Place the child on the scale, which will then display the measured weight.
4. The measured value on the display flashes 3 times and then remains displayed.
5. Turn off the scale by holding down the **ON/OFF/TARE** button, or the scale will automatically turn off after 2 minutes.

FIRST USE, WEIGHING WITH THE APP

- The scale can also be used in conjunction with a smartphone app to help you monitor your child's development and growth using measurement history and graphs.

1) Download the app and register

- a. Scan the QR code or enter the URL link and install the app.
- b. Fill in your personal information and register.

truellife.eu/58605/app

2) Pair the scale with the app

- a. Turn on the app and activate Bluetooth on your smartphone.
- b. In the app, click the **Add device** button or the + in the upper right corner.
- c. From the menu select **Baby scale (Bluetooth)**.

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